Operating Instructions

Please read carefully and keep this copy for further use.

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OSSEOSTAP INSTRUMENTS (EN)

Operating Instructions and safety recommendations for the use of OSSEOSTAP rotary instruments in the medical field.

Bien-Air Surgery SA rotary instruments are made of high-grade materials. They are designed, constructed and produced with the utmost care. These quality instruments assure best working results and long service life provided they are used in the proper manner. Therefore, the following instructions for use and safety recommendations have to be observed.

Improper use can lead to damage to the tissue, premature wear, destruction of the instruments and injury to the operator, patient or other persons.

OSSEOSTAP rotary instruments are supplied « sterile ».

Do not use an instrument if its packaging is open, damaged or after expiry date indicated on the packaging because there is no more protection against the risk of contamination.

1. Proper use

- Observe this information and the instructions for use for Bien-Air Surgery SA OSSEOSTAP system when using these rotary instruments.
- Prior to use, make sure that the instruments are in perfect technical condition and that sterility is assured. Never use damaged instruments.
- Insert rotary instruments in the hand piece until the stop position is reached. (see OSSEOSTAP operating instructions)
- Prior to operation, make sure that the instrument is firmly seated.
- Test each instrument at the desired speed before using it. If it vibrates or functions in an abnormal manner, check that the instrument is installed properly or replace it.
Never use an instrument which is not operating correctly because it might otherwise be the cause of cuts or wounds.
- Instruments diameter color coding:
In order to allow an easy and fast identification of the various OSSEOSTAP instrument diameters, a colored ring is affixed on the shank. The correspondences are the following:

Violet	Ø 0.5 mm
Blue	Ø 0.6 mm
Green	Ø 0.7 mm
Yellow	Ø 0.8 mm
Orange	Ø 1.0 mm
Grey	Ø 1.4 mm
Brown	Ø 1.8 mm
Red	Ø 2.3 mm
Black	various Ø / others than above

Bien-Air Surgery SA provides a small plate in order to preserve the color coding of the diameters near tools (ref. 1304358 (5 pieces per box))

2. Intended use

The OSSEOSTAP instruments are intended for the light drilling of ossicle bones.

3. Indication for use

Product intended for professional use. The OSSEOSTAP instruments have been designed for the light drilling of bones as part of surgical operations such as stapedotomy or ossiculoplasty.
The selection of the correct instrument, especially the correct drilling diameter for stapedotomy prostheses, is the responsibility of the surgeon.

4. Contraindications

None currently known.

5. Operation

Due to its slim geometry allowing access to the middle ear, the OSSEOSTAP handpiece is not adapted for important drillings. Excessive contact pressure has to be avoided.

- It can lead to damage the instrument and thus generating an important heating of the bur guide
- At the worst, breakage of the instrument can occur if too high contact pressure is applied. Jamming or leverage of the instrument must be avoided during operation.

6. Cleaning, disinfection

After each use, the instruments must be pre-cleaned, cleaned, inspected and sterilized.

Point of use cleaning


Immediately, after the surgical act, proceed as follows:

- Remove the instrument from the handpiece.
- Remove blood, tissue or bone residues (risk corrosion).
- Never submerge the OSSEOSTAP instruments in physiological salt water solutions (NaCl solution), because corrosion risks.
- Ensure that the OSSEOSTAP instruments do not dry before manual or automatic cleaning by wrapping them in non-woven towelettes (pre-soaked in water)

- Manual or automatic cleaning must be done max 4h after the point of use.

Manual cleaning

- Rinse the instrument under running tap water (cold). While rinsing, move the distal bushing back and forth.
- Brush the tip and both bushings thoroughly, parallel to the cutting flutes, with a nylon soft bristles for at least 10 seconds. Continue brushing until soil is no longer visible.

- Soak the instrument in ultrasonic bath with enzymatic pH neutral detergent* (like Steris Prolystica® 2X Concentrate Enzymatic Presoak and Cleaner) solution for at least 3 minutes.
- Remove instrument from detergent bath and rinse thoroughly under running tap water (cold) for at least 10 seconds.
- Blow dry with pressurized air.

* Time/temperature/concentration according to detergent manufacturer's instructions.

Or

Automatic cleaning

- Rinse the instrument under running tap water (cold).
- Brush the instrument with nylon soft bristles for at least 10 seconds or until soil is no longer visible.
- Place the instruments in the appropriate washer/disinfectant basket and treat via a standard instrument washer/disinfectant cycle.
- Exclusively use a washer/disinfectant validated as per (ISO 15883).

This product has been validated with an automated cycle consisting of the following stages:

1. Pre-wash: Cold tap water (<45°C / 113°F) for minimum 2 minutes.
2. Washing: Hot tap water, 50°C to 60°C (122°F to 140°F) with enzymatic pH neutral detergent (like Steris Prolystica® 2X Concentrate Enzymatic Presoak and Cleaner), minimum 5 minutes.
3. Neutralization: Cold tap water (<45°C / 113°F) for minimum 2 minutes.
4. Rinsing: Cold critical water acc to AAMI TIR34 (<45°C / 113°F) for minimum 2 minutes.
5. Thermal disinfection by rinsing: Hot critical water acc to AAMI TIR34 (90°C / 194°F) for minimum 5 minutes. Operator is responsible for the implemented value A0 according A0 concept described in EN ISO 15883 (For example, A0 600 90°C (194°F) / 1min.).
6. Dynamic drying under ventilation: at 70°C (158°F), for minimum 22 minutes.

Comments:

- Comply with the washer/disinfectant loading instructions provided by the manufacturer.
- Make sure that all the instruments have been correctly attached to the baskets. Make sure that the instruments do not touch one another.
- Remove the instruments from the washer or disinfectant immediately after the machine stops to avoid corrosion.
- It is recommended to use anticorrosive disinfection and cleaning agents. Contact with H₂O₂ (hydrogen peroxide), chlorox disinfection and cleaning agents as well as oxalic acids are to be avoided because these agents may lead to pitting and corrosion.

7. Cleaning check, sorting out of dull instruments and lubrication

After each cleaning, carefully inspect the instrument with a magnifying glass to make sure that all visible contamination has been eliminated. If contamination is present, repeat the cleaning process.

The lifespan of the OSSEOSTAP instruments is basically determined by the wear and damages due to the use.

After each cleaning, and before sterilization, inspect the instrument to make sure that it is not dull or damaged.

Look for the following indications:

- dull (rounded) blades and breakouts
- blank spots on diamond instruments
- damaged or bent shank
- bushings without trace of wear, one must be mobile axially, the second must be fixed

Under no circumstances use dull or damaged instruments any more. In the event of light bending of the shaft, straighten it by delicately drawing it between 2 fingers.

Non observance of this procedure may generate instrument breakage and risk of injury to the patient and operator.

Following inspection and before sterilization, lubricate the instrument by applying a light coating of medical oil (such as Bien-Air "Lubrifiant", REF 1600064) on the bushing area. Move the mobile bushing forward and backward at least 2 times.

8. Sterilization

Sterilization by steam/moist heat is recommended. The following sterilization parameters, using a Pre-Vac cycle, are recommended by Bien-Air Surgery SA to provide a sterility assurance level (SAL) of 10⁻⁶ of wrapped instruments (In the USA, use an FDA approved sterilization wrap).

Temperature	132°C (270°F)	134°C ¹ (273°F)	134°C ¹ (273°F)	135°C (275°F)
Time	4 min.	3 min.	18 min. ²	3 min.
Drying time	30 min. ³			

¹ Not for users in US healthcare facilities.

² Parameters recommended by the World Health Organization for treating instruments in the event of contamination by Non-Conventional Transmissible Agents (NCTA). Parameters recommended by Bien-Air Surgery SA.

³ Refer to sterilizer manufacturer recommendations for drying times per load configuration.

After sterilization, let the devices cool down to room temperature without forced cooling.

Do not exceed a temperature of 138°C (280°F) the guidance bushings will get damaged.

9. Storage / Environmental conditions

The instruments are to be stored in an environment free of dust and humidity guaranteeing no recontamination.

Work:

- Temperature +10°C to +30°C (+50°F to 86°F)
- Relative humidity 20% to 80%, including condensation
- Atmospheric pressure 700 hPa to 1060 hPa.

Transport:

- Temperature -25°C to +70°C (-13°F to 158°F)
- Relative humidity 10% to 100%, including condensation
- Atmospheric pressure 600 hPa to 1060 hPa.

Storage:

- Temperature +10°C to +30°C (+50°F to 86°F)
- Relative humidity 20% to 80%, including condensation
- Atmospheric pressure 600 hPa to 1060 hPa.

10. Safety and liability

- The user is responsible to check the product prior to use whether it is suitable for the intended purpose.
- The user is responsible for the application of the instruments.
- In case of contributory negligence by the user, Bien-Air Surgery SA partially or totally declines liability for all resulting damages, particularly if these are due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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