


E.N.T.

RIGID AND FLEXIBLE ENDOSCOPES

ENDOSCOPES

SINUSCOPES / LARYNGOSCOPES / OTOSCOPE



E.N.T.


RIGID AND FLEXIBLE ENDOSCOPES

We offer exclusive design, excellent quality and cost effective endoscopes.
 Outstanding image quality. • Exceptional luminosity and brightness.
 Storz / Wolf and ACMI connections.




SINUSCOPES

DIAMETER	3.0mm		
WORKING LENGTH	Approx. 140mm		
DIRECTION OF VIEW	0°	30°	70°




SINUSCOPES

DIAMETER	4.0mm			
WORKING LENGTH	Approx. 175 mm			
DIRECTION OF VIEW	0°	30°	45°	70°



OTOSCOPES

DIAMETER	3.0mm			
WORKING LENGTH	100 mm	140 mm		
DIRECTION OF VIEW	0°	30°	0°	30°



LARYNGOSCOPES

DIAMETER	6.0mm	8.0mm
WORKING LENGTH	Approx. 175 mm	Approx. 175 mm
DIRECTION OF VIEW	70°	70°



FLEXIBLE SCOPE

FIELD OF VIEW	85°	
DEPTH OF FIELD	2,5 - 50 mm	
OUTER DIAMETER INSERTION TUBE	3.4 mm	
OUTER DIAMETER DISTAL END	3.4 mm	
WORKING LENGTH	300 mm	
BENDING RANGE OF DISTAL TIP	up 130°	down 130°

Rev.06.11.18

ECLERIS RESERVES THE RIGHT TO MODIFY THIS PRODUCT WITHOUT PRIOR NOTICE