

PROCAM SD

VIDEO CAMERA

- ENDOSCOPY / COLPOSCOPY / MICROSCOPY



- High resolution one CCD camera
- Less than one lux sensitivity
- Ideal for the office, affordable
- Great image quality
- Optimize light control
- Ergonomic head with 3 controls:
 - White balance
 - Windows function (two levels)
 - Anti-moiré function (two levels)

TECHNICAL SPECIFICATIONS

Resolution	Pal 625 lines / ntsc 525 lines
Ccd	1/3"
Shutter speed	Up to 1/100.000 sec.
Min. Illumination	2 lux f1.4
Video output	Composite / S-VIDEO
Power supply	100v / 240v
Frequency	50 - 60hz