



EVERLUX UNIVERSAL

Fuente de Luz LED súper brillante con conexión de fibra óptica universal

LED 80W. 50,000 horas de vida.
Luz fría, blanca, de alta intensidad.
Dimmer electrónico.
Conexión de fibra óptica universal.

TECHNICAL SPECIFICATIONS

LED	80 W
Light color temperature	6000° k
Power consumption	175 watts
Class	I
Type	BF
LED life	50.000 hours approx.
Over temp. protection	Automatic LED system power-off
Dimmer	Electronic
Fiber optic connector	Universal, compatible with major brands
Dimensions	202 X 290 X 95mm (width x depth x height)
Weight	2.5 Kg.
Power supply	100- 240 V
Frequency	50 / 60Hz

EVERLUX

Fuente de luz LED súper brillante



Compacta y portátil.
LED 80W. 50,000 horas de vida.
Luz fría, blanca, de alta intensidad.
Conexión individual para la mayoría de los cables OF (Storz, Olympus, tipo Wolf).

TECHNICAL SPECIFICATIONS

LED	80 W
Light color temperature	6000° k
Power consumption	175 watts
Class	I
Type	BF
LED life	50.000 hours approx.
Over temp. protection	Automatic LED system power-off
Dimmer	Electronic
Fiber optic connector	Individual, compatible with major brands
Dimensions	202 X 290 X 95mm (width x depth x height)
Weight	2.5 Kg.
Power supply	100- 240 V
Frequency	50 / 60Hz

EVERLUX PRO

Fuente de luz LED súper brillante



Luz fría, blanca, de alta intensidad.
LED 80W. 50,000 horas de vida.
Conector integrado para cámara.
Conexión de fibra óptica compatible con las principales marcas.

TECHNICAL SPECIFICATIONS

SERIES	LED LIGHT SOURCE
LED	80 W
LIGHT COLOR TEMPERATURE	5900°K
POWER SUPPLY VOLTAGEN*	90 - 264 V
OF PHASES	1 (SINGLE-PHASE)
FREQUENCY	50 / 60Hz
POWER CONSUMPTION	175 W
CLASS	I
TYPE	B
LED SYSTEM SERVICE LIFE	50.000 HOURS APPROX.
REFRIGERATION	FORCED THROUGH A FAN
OVER TEMP. PROTECTION	AUTOMATIC LED SYSTEM POWER-OFF
DIMMER	ELECTRONIC
FIBER OPTIC CONNECTOR	BUILT-IN REVOLVING TURRET WITH MULTIPLE ADAPTERS
DIMENSIONS	202 X 290 X 95mm (width x depth x height)
WEIGHT	3.0 Kg.

Accesorios: Cable de fibra óptica

Largo: 1,80 metros
Conector: Storz-Storz

