

# HD SPLITTER

Full 1920 x 1080 p.



- In-line beam splitter for superior balance.
- Ergonomic and compact design.
- Compatible with Microscopes and Colposcopes of multiple brands.
- Outstanding image quality.

# HD SPLITTER

Full 1920 x 1080 p.

**This beam splitter integrates a full HD video camera in an ergonomic and compact way.**

The device is installed between the optical head and the binocular. This configuration results in great balance since it avoids the use of photographic and/or video equipment located on the side of the optical head which usually alter the stability of the movements of it.



This product is complemented by sophisticated image capture systems from the ECLERIS ENDODIGI family.

ECLERIS RESERVES THE RIGHT TO MODIFY THIS PRODUCT WITHOUT PRIOR NOTICE