

STROBOLIGHT XENON

- Both continuous and stroboscopic light with same xenon bulb
 - The brightest stroboscope in the market
- Automatic voice sincronization • Total footswitch control



STROBOLIGHT XENON



- Our Strobolight stroboscope is developed with an innovative electronic light control system (Xenon Auto Transition Technology) that provides exceptional light flow. It has the unique feature in the market that allows functioning with the same 180-watts Xenon bulb both in continuous and stroboscopic mode; this feature makes it extremely comfortable and simple to use.
- Contact microphone, automatic voice detection and automatic change of illumination mode, from continuous to stroboscopic.
- Full footswitch operation.
- Automatic detection and synchronization with voice frequency.
- Provided with manual function for pathological patients.
- Includes a rotary turret that fits all main fiber optic connections available on the market (Storz, ACMI, Wolf and Olympus).
- Amplified audio output.

ECLERIS RESERVES THE RIGHT TO MODIFY THIS PRODUCT WITHOUT PRIOR NOTICE

Rev. 16.01.15