

RHYNOSOFT

COMPUTARIZED RHINOMANOMETER

FAST AND ACCURATE



RHYNOSOFT

COMPUTARIZED RHINOMANOMETER

- Customized test
- Easy to install
- Patient database
- Automatic calculation of main parameters



TECHNICAL SPECIFICATIONS

CHANNELS SENSIBILITY

2 independent (pressure and flow)
Pressure channel: 3Pa (0,3mm H₂O)
Flow channel: 5cm³/S

APPLICATIONS

Automatic calculation of the air flow resistance.
Automatic calculation of respiratory
Automatic calculation of pressure and flow inhalation
and exhalation peaks
Automatic flow / pressure graphic generation

POWER SUPPLY PC CONNECTION

110V / 220V
USB port

MINIMAL REQUIREMENTS OPERATING SYSTEM

Pentium II 400 MHz, 64MB RAM, 1G of free disk space
Windows®98SE, Windows® XP, Windows® ME

ECLERIS RESERVES THE RIGHT TO MODIFY THIS PRODUCT WITHOUT PRIOR NOTICE